



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Product Manual of DAQIDA

Choking Rescue Device set – Dual-Tube Manual Suction System for Adults and Children

B2B Product Specification & Description Manual | DAQIDA

Product Overview

The Choking Rescue Device Set – Dual-Tube Manual Suction System for Adults and Children is an OEM-ready, FBAO emergency clearance kit designed for wholesale and global B2B distribution. As a vital second-line airway clearance device, this portable anti-choking suction pump requires no electricity or batteries, making it an essential wholesale safety product for medical distributors, school safety programs, and household emergency kits.



Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Key Product Features

Seal Cover: Tightens to ensure the internal plunger system remains hermetically sealed for optimal compression.

Handle & Plunger: Ergonomically designed; pulled back to create powerful negative pressure, efficiently extracting airway obstructions.

Holder & Vent Valve: Press and hold to establish an airtight facial seal, expelling internal air while preventing backflow to guarantee maximum suction power.

Face Mask: Ultra-soft, skin-friendly medical PVC material designed to fit facial contours comfortably without risking skeletal trauma. Includes interchangeable Adult Mask and Pediatric Mask.

Open Tube: An optional auxiliary tool used to keep the airway or mouth accessible when required during complex rescue scenarios.



How to Use

Step 1: PLACE: Cover the patient's mouth and nose with the PVC mask, pressing firmly to establish an airtight facial seal.

Website: www.daqida.com **WhatsApp / WeChat:** +86 187 5767 2042

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Step 2: PUSH: Grip the PP+TPE dual-tube body and push downward to discharge internal air safely through the one-way valve.

Step 3: PULL: Keep the mask sealed and pull upward to generate powerful negative pressure that extracts the foreign object.

Step 4 REPEAT: If needed, maintain the seal and repeat the push-pull cycle until the airway is

EASY TO USE IN CHOKING EMERGENCIES

1 PLACE

2 PRESS

3 PULL

Place the mask over the victim's nose and mouth.

Press and hold the holder&vent valve with your hand so that it forms a good all-around seal on the victim's face. This allows it to generate strong suction.

Pull back the handle&plunger while securing the mask to the victim's face. This will create suction. If the obstruction does not come out on the first pull, continue pushing and pulling the handle&plunger without removing the mask from the victim's face until the obstruction is removed.

▼ PLACE



▼ PRESS



▼ PULL



FOR OTHERS



FOR YOURSELF





This product uses mechanical action to create one-way suction and is designed with a one-way valve to prevent additional injury. It's made of high-quality material, is portable and durable, reusable, and doesn't need replacement parts. It's not intended for medical use, but it can be used in emergency situations and is easy to use with clear instructions.

cleared.

Materials & Technical Selection & Engineering

Main Body (PP + TPE): Combines rigid Polypropylene (PP) for structural stability with flexible Thermoplastic Elastomer (TPE) for smooth, fatigue-free operation.

Face Mask (Medical-Grade PVC): Biocompatible, transparent, and soft medical-grade PVC, allowing real-time visual monitoring during emergency response.

Safety Mechanism: Purely mechanical one-way suction action that safely expels foreign objects without pushing them deeper or requiring batteries.

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

MORE CUSTOMIZATION

**PLEASE CLICK ON
THE BOTTOM OF THE SCREEN**



QUALITY
ASSURANCE



OEM/ODM
SUPPORT



FAST
DELIVERY



PROFESSIONAL
SERVICE



Click this button to immediately
get factory prices and more services

Application Scenarios & Use Cases

Household Emergency Readiness: An essential safety device for family living rooms, kitchens, and bedrooms to protect household members of all ages from unexpected choking emergencies.

Automotive & Travel Safety: Compact and portable design ideal for storage in vehicle glove boxes, travel bags, or road trip emergency kits for on-the-go protection.

School & Institutional Protection: A reliable, battery-free safety addition for classrooms, cafeterias, daycares, and campus health stations where quick, non-invasive airway clearance is critical.

Public Dining & Hospitality: Recommended for restaurants, hotel dining halls, and hospitality venues to support staff response readiness during dining incidents.

Professional & Commercial Distribution: Perfect for medical distributors, safety equipment wholesalers, and e-commerce brands seeking high-demand, OEM-ready household safety gear.

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

READY WHEN YOU NEED IT



Technical Specifications

Product Specifications

Product Name	Choking Rescue Device
Set Specification	Single Pack
Product Color	Transparent Yellow, Transparent Black, Transparent Blue, Main Body Yellow
Product Size (cm)	Adult Mask: 12.5×10.0×7.5CM, Kids Mask: 10×7.5×6CM, Main Body: 16×9.5×9.5CM
Product Package Size (cm)	30×13×7CM
Product Material	PE+PVC
Product Net Weight (g)	165.8 g
Product Gross Weight (with Package, g)	246.5 g
Packaging Material	Yellow Bag Version
Outer Carton Size (cm)	52×48×52CM
Net Weight per Carton (kg)	12.3 KG
Gross Weight per Carton (kg)	—
Quantity per Carton (Sets)	50 PCS
Custom Color MOQ (Sets)	—
Custom Logo MOQ (Sets)	—
Invoicing	(Please inform in advance) 1% Regular Invoice / 13% VAT Invoice

Website: www.daqida.com **WhatsApp / WeChat:** +86 187 5767 2042

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Patented Choking Rescue Device Set – Dual-Tube Manual Suction System

Our production lines feature official United States Design Patents and EUIPO registrations for our proprietary structural designs.

Sourcing a patented dual-tube choking rescue device protects your supply chain from market infringement and ensures full legal compliance.

As an experienced OEM anti-choking device manufacturer, we combine independent patents with FDA compliance to deliver reliable emergency solutions.



Website: www.daqida.com WhatsApp / WeChat: +86 187 5767 2042

Email: sean@daqida.com Address: Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province