



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Product Manual of DAQIDA

Motion sickness glasses

B2B Product Specification & Description Manual | DAQIDA

Product Overview

As a direct B2B wholesale supplier, DAQIDA produces high-performance anti car sick glasses engineered to prevent motion discomfort during travel. Commonly used as car sickness glasses and seasick glasses



Key Product Features

Effective Motion Relief: Rapidly mitigates nausea, dizziness, and discomfort caused by cars, ships, planes, or trains within 10–12 minutes of wear.

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Drug-Free & Reusable: A 100% natural alternative to motion sickness pills, featuring no side effects or drowsiness for adults and children.

Universal Compact Fit: Lightweight, lensless structure fits comfortably over most facial shapes or prescription eyewear.

ANTI-MOTION SICKNESS GLASSES
Relieve motion sickness during travel

Swayable Liquid
Blue liquid, making the eyes feel objects moving in multiple directions, Reducing the nerve difference where the ears perceive object movement while the eyes perceive object stillness.

Lightweight, no pressure.
Made of imported PC material, the whole machine weighs only 53g.

53g

Snap design for easy installation and removal.
Four glasses frames around all feature snap fastener design for easy installation and removal.

Swayable Liquid | **Lightweight Comfort** | **Safe & Reliable** | **Easy to Use**

How to Use

Put On at First Sign of Discomfort: Wear the glasses immediately when motion sickness or nausea begins during travel.

Wear for 10 to 12 Minutes: Look at a stationary object (like a phone or book); the liquid rings will adjust to restore visual equilibrium.

Remove Once Relieved: Take off the glasses once symptoms subside and enjoy a comfortable journey.

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Materials & Technical Selection & Engineering

High-Impact ABS Frame: Molded from durable, lightweight ABS resin to prevent cracking or deformation during outdoor travel.

Skin-Safe Liquid Silicone: Utilizes soft, hypoallergenic silicone for nose pads and temple arms, ensuring a pressure-free fit without skin irritation.

Chamber Liquid Technology: Features sealed circular rings filled with colored liquid that move dynamically to simulate a horizon.



Application Scenarios & Use Cases

Automotive & Road Trips: Primary car sickness glasses for passengers prone to nausea while reading or using mobile devices in cars.

Maritime & Cruise Travel: Essential seasick glasses for boat rides, fishing trips, and ocean cruises.

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Retail & Travel Accessories Supply: Ideal inventory for airport travel shops, pharmacies, and Amazon e-commerce brands.

Wholesale & Private Labeling: Factory-direct motion sickness glasses available for B2B buyers seeking custom logo branding and bulk packaging.



Technical Specifications

Product Name	Anti Motion Sickness Glasses
Package Specification	Single Pack
Product Color	White
Product Size (CM)	16×6×4cm
Product Package Size (CM)	17.6×8.4×4.5cm
Product Material	Plastic
Product Weight (Net)	42.9g

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province



DAQIDA | Yiwu Daqida Plastic Products Co., Ltd.

Product Weight with Package (Gross)	75.7g
Packaging Material	Color Box
Outer Carton Size (CM)	52×48×52cm
Carton Net Weight (KG)	11.4kg
Carton Gross Weight (KG)	13kg
Quantity per Carton (Sets)	150 sets
MOQ for Custom Color (Sets)	TBD
MOQ for Custom Logo (Sets)	TBD
Invoicing	1% ordinary VAT invoice / 13% VAT invoice (please inform in advance)

Website: www.daqida.com **WhatsApp / WeChat:** [+86 187 5767 2042](https://wa.me/8618757672042)

Email: sean@daqida.com **Address:** Building 1, No. 47 Longxiang Road, Nansanli, Yiwu City, Zhejiang Province